
## Revisions

SO 0191-5


Rev from	Rev to	ECO #
1204	2010-01	5253-10

### Notes:


- BD Cat. Number 220231, 220232, 220233, 220234
- Blank (Sheet) Size : Length: 23 5/8" Width: 3 15/32"  
Number of Pages: 2 Number of Sheets: 1  
Page Size: Length 23 5/8" Width 3 15/32" Final Folded Size: Special Fold\*
- Style (see illustrations below): # 1


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number OEM - Copan for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No  
No. of Colors: 1 PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

\*Peel-off insert / Label. Insert perforated 3/8" from bottom with 3/16" notch for easy peel-off. Top insert size 4 3/16" x 3 15/32" includes 3/8" bottom portion of glued area. Insert is folded under with one horizontal fold and followed by two parallel folds to form the size of 3 1/4" x 3 15/32." Insert to be glued on top 5/16" and bottom 3/8" to roll label, size 3 15/32" (web) x 4 3/8".

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 <b>BD</b> Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA
Proofer	Date		
Checked By	Date		
Part Number:  L009329/H283	Category and Description  Peel-off Insert / Label, CultureSwab MaxV	Sheet: 1 of 3	A
		Scale: N/A	

Pokyny vám poskytnú miestni zástupcov spoločnosti BD. / Kontakt den lokale BD representant for at få instruktioner. / Kasutusi juhiste sahtes kontakteeruge oma kohaliku BD esindajaga. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / A használati utasítást kérje a BD helyi képviselőjétől. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcję użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o seu representante local da BD para obter instruções. / Instrukcie zyskate u miestneho zástupcu spoločnosti BD. / Kontakta lokal Becton Dickinson-representant för anvisningar. / Свържете се с местния представител на BD за инструкции. / Contactați reprezentantul dumneavoastră local BD pentru instrucțiuni. / Talimatlar için yerel BD temsilcilerinize danışın. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Для получения инструкций свяжитесь с местным представителем компании BD. / Өзіндік жергілікті БД өкіліне жүгініп нұсқау алыңыз. / Kontaktiraj lokalnog predstavnika BD za upute.

**INTENDED USE**

**BD CultureSwab™** MaxV devices are sterile ready-to-use systems intended for the collection, transport and preservation of clinical specimens for bacteriological examination.

**SUMMARY AND PRINCIPLES**

One of the routine procedures in the diagnosis of bacterial infections involves the collection and safe transportation of a clinical specimen from the patient to the laboratory. This can be accomplished using the **BD CultureSwab** MaxV collection and transport device. Each **BD CultureSwab** MaxV unit is comprised of a sterile peel pouch containing a rayon-tipped swab applicator used to collect the sample and a tube containing transport medium into which the swab applicator is placed after sampling.

The **BD CultureSwab** MaxV transport media (Amies Liquid Medium and Liquid Stuart Medium) are non-nutritious, buffered with phosphate and provide a reduced environment due to their formulation with sodium thioglycollate.<sup>1</sup> Organisms in the sample material are protected from drying by moisture in the transport medium. The medium is designed to maintain the viability of organisms during transit to the laboratory. **BD CultureSwab** MaxV pouches are made of a plastic film which retards the penetration of atmospheric air into the product.

**BD CultureSwab** MaxV media are contained in a constricted (venturi) tube. Nitrogen gas is flushed into the transport tube during the media filling and capping process. During final packaging of the swab and tube, air is removed from the pouch by vacuum and nitrogen gas is flushed inside.

**REAGENTS**

The nominal formula for each medium per L of distilled water is:

Amies Liquid Medium		Stuart Liquid Medium	
Calcium Chloride	0.10 g	Calcium Chloride	0.10 g
Disodium Phosphate	1.15 g	Mercaptoacetic Acid	1.0 mL
Magnesium Chloride	0.10 g	Sodium Glycerophosphate	10.00 g
Monopotassium Phosphate	0.20 g		
Potassium Chloride	0.20 g		
Sodium Chloride	3.00 g		
Sodium Thioglycollate	1.00 g		

**Precautions:** For *in vitro* Diagnostic Use.

It must be assumed that all specimens contain infectious microorganisms; therefore, all specimens should be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste.

ⓧ **BD CultureSwab** MaxV is for single use only; reuse may cause a risk of infection and/or inaccurate results.

**Storage:** Store **BD CultureSwab** MaxV devices at 5–25°C.

**Product Deterioration:** Contents sterile if unopened or not damaged. Do not use if they show evidence of damage, dehydration or contamination. Do not use if past expiration date.

**SPECIMEN COLLECTION AND HANDLING**

For specific recommendations about collection of specimens for microbiological analysis and primary isolation techniques, consult appropriate references.<sup>2-4</sup>

Once a swab sample is collected, it should be placed in the tube of medium, transported to the laboratory as soon as possible and cultured onto appropriate primary isolation media.

**PROCEDURE**

**Materials Provided:** Fifty (50) units of sterile **BD CultureSwab** MaxV devices, which contain media, are contained in each Vi-Pak Pouch.

**Materials Required But Not Provided:** Appropriate materials for isolating, differentiating and culturing aerobic and anaerobic bacteria. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations.

**Directions For Use:**

The directions for use are printed on each **BD CultureSwab** MaxV unit, along with descriptive diagrams, and are summarized as follows:

1. Peel open the **BD CultureSwab** MaxV pouch.
2. Remove cap from transport tube.
3. Remove applicator swab and collect specimen.  
During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.
4. Place applicator swab in transport tube.
5. Record patient's name and information on tube label.
6. Send specimen to the laboratory for immediate analysis.

**EXPECTED RESULTS**

The survival of bacteria in a transport medium depends on many factors. These include the types of bacteria, duration of transport, storage temperature, concentration of bacteria in the sample and formulation of the transport medium. **BD CultureSwab** MaxV devices maintain viability of many microorganisms for 24–48 h. For fastidious bacteria, such as *Neisseria gonorrhoeae* and *Streptococcus pneumoniae*, swab specimens should be plated directly onto culture medium or transported immediately to the laboratory and cultured within 24 h.

**LIMITATIONS OF THE PROCEDURE**

**BD CultureSwab** MaxV Amies Liquid Medium and Stuart Liquid Medium are intended for the collection and transport of bacteriological samples only. Preferred samples for anaerobic investigations are: tissue samples obtained during surgical procedures, biopsies from tissue or bone, fluid, pus, or aspirates collected using a syringe. For detailed information and recommendations for transporting fluid and tissue specimens for anaerobic culture, refer to specific publications.<sup>3-7</sup> Samples containing viruses or chlamydia should be collected and transported using alternative specific transport systems.

**BD CultureSwab** MaxV has not been validated with direct antigen tests. For approved collection methods for direct antigen tests, refer to the appropriate manufacturers' kit insert.

**PERFORMANCE CHARACTERISTICS**

Recovery studies were performed using **BD CultureSwab** MaxV products with a variety of aerobic organisms. Swabs were dosed with inoculum and inserted into the transport tube containing transport media. The tubes were stored at room temperature prior to subculturing onto appropriate culture media. Listed are the organisms tested with the various transport media.

Organism	Transport Media	
	Amies Liquid form	Stuart Liquid form
<i>Neisseria gonorrhoeae</i> ATCC™ 43069	*	*
<i>Streptococcus pyogenes</i> ATCC 19615	*	*
<i>Haemophilus influenzae</i> ATCC 10211	*	*
<i>Streptococcus pneumoniae</i> ATCC 6305	*	*

**AVAILABILITY** - See end of insert.

**REFERENCES**

1. Amies CR. 1967. A modified formula for the preparation of Stuart's transport medium. Can. J. Public Health. 58: 296-300.
2. Isenberg H.D., F.D. Schoenkencht, and A. von Graevenitz. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating ed., S.J. Rubin. American Society for Microbiology, Washington, D.C.
3. Balows, A., W.J. Hausler, Jr, K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy (ed.). 1991. Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
4. Isenberg, H.D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. Zavala, M.K., D.M. Citron, E.J.C. Goldstein. 1998. Evaluation of a novel specimen transport system for anaerobic bacteria. Clin. Infect. Dis. 25 (supplement 2): S132-S133.
6. Perry, J.L. 1997. Assessment of swab transport systems for aerobic and anaerobic organism recovery. J. Clin. Microbiol. 35: 1269-1271.
7. Summanen, P., E.J. Baron, D.M. Citron, C.A. Strong, H.M. Wexler, and S.M. Finegold. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Co., Belmont, Calif.

**TECHNICAL INFORMATION:** In the United States, telephone Technical Services, toll free (800) 638-8663.

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo and CultureSwab are trademarks of Becton, Dickinson and Company. © 2010 BD.

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